

For Research Use Only

FITC Plus Anti-Mouse CD41 (MWReg30)

Catalog Number: FITC-65189

1 Publications



Basic Information

Catalog Number:

FITC-65189

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120493

GeneID (NCBI):

16399

UNIPROT ID:

Q9QUM0

Full Name:

integrin alpha 2b

Purification Method:

Affinity purification

CloneNo.:

MWReg30

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Longxing Miao	39172150	Naunyn Schmiedebergs Arch Pharmacol	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse peripheral blood platelets were surface stained with 0.5 ug FITC Plus Anti-Mouse CD41 (MWReg30) (FITC-65189, Clone: MWReg30) (red) or FITC Plus Rat IgG1 Isotype Control (HRPN) (FITC-65212, Clone: HRPN) (blue). Cells were not fixed.