

For Research Use Only

# FITC Plus Anti-Mouse CD4 (RM4-4)

Catalog Number: FITC-65199

1 Publications



## Basic Information

<b>Catalog Number:</b> FITC-65199	<b>GenBank Accession Number:</b> BC039137	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100ug , 500 µg/ml	<b>GeneID (NCBI):</b> 12504	<b>CloneNo.:</b> RM4-4
<b>Source:</b> Rat	<b>Full Name:</b> CD4 antigen	<b>Excitation/Emission maxima wavelengths:</b> 495 nm / 524 nm
<b>Isotype:</b> IgG2b, kappa		

## Applications

**Tested Applications:**  
FC

**Cited Applications:**  
FC

**Species Specificity:**  
Mouse

**Cited Species:**  
mouse

## Background Information

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

## Notable Publications

Author	Pubmed ID	Journal	Application
Wenlong Lu	37370477	Animals (Basel)	FC

## Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with CoraLite® Plus 647 Anti-Mouse CD3 (CL647-65077, Clone: 17A2) and 1 ug FITC Plus Anti-Mouse CD4 (FITC-65199, Clone: RM4-4) or FITC Plus rat IgG2b isotype control. Cells were not fixed.