

For Research Use Only

FITC Rat IgG2a Isotype Control (2A3)



Catalog Number: FITC-65209

Basic Information

Catalog Number:

FITC-65209

Size:

100ug, 0.5 mg/mL

Source:

Rat

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Affinity purification

CloneNo.:

2A3

Applications

Tested Applications:

FC

Species Specificity:

Background Information

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.25 ug FITC Rat IgG2a Isotype Control (FITC-65209, Clone: 2A3) (blue) or 0.25 ug FITC Anti-Mouse CD4 (FITC-65141, Clone: RM4-5) (red). Cells were not fixed.