

For Research Use Only

FITC Plus Anti-Rat CD43 (W3/13)

Catalog Number: FITC-65237



Basic Information

Catalog Number: FITC-65237	GenBank Accession Number: GeneID (NCBI): 24796	Purification Method: Affinity purification
Size: 100ug, 500 µg/ml	Full Name: sialophorin	CloneNo.: W3/13
Source: Mouse	Calculated MW: 38 kDa	Excitation/Emission maxima wavelengths: 495 nm / 524 nm
Isotype: IgG1, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Rat

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 Wistar rat splenocytes were surface stained with 0.5 μ g FITC Plus Anti-Rat CD43 (FITC-65237, Clone: W3/13) or 0.5 μ g Control Antibody. Cells were not fixed.