

For Research Use Only

FITC Plus Anti-Mouse CD49b (HMA2)

Catalog Number: FITC-65245



Basic Information

Catalog Number:
FITC-65245

Size:
100ug, 500 µg/ml

Source:
Armenian Hamster

Isotype:
IgG

GenBank Accession Number:
BC065139

GeneID (NCBI):
16398

Full Name:
integrin alpha 2

Purification Method:
Affinity purification

CloneNo.:
HMA2

Excitation/Emission maxima
wavelengths:
495 nm / 524 nm

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with APC Anti-Mouse NK1.1 (CD161) (APC-65138, Clone: PK136) and 0.5 ug FITC Plus Anti-Mouse CD49b (FITC-65245, Clone: HMA 2) or 0.5 ug FITC Plus Armenian Hamster IgG. Cells were not fixed.