

For Research Use Only

FITC Plus Anti-Mouse GL7 (GL-7)

Catalog Number: FITC-65261



Basic Information

Catalog Number:

FITC-65261

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgM, kappa

GenBank Accession Number:

GeneID (NCBI):

104230

Full Name:

lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



1X10⁶ unstimulated or anti-CD3/CD28 stimulated mouse splenocytes were surface co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) and 0.5 ug FITC Plus Anti-Mouse GL7 (FITC-65261, Clone:GL-7) or 0.5 ug Isotype Control. Cells were not fixed.