

For Research Use Only

# FITC Plus Anti-Rat CD43 (W3/13) Rabbit Recombinant Antibody



Catalog Number: FITC-65556

## Basic Information

Catalog Number:  
FITC-65556  
Size:  
100ug , 500 µg/ml  
Source:  
Rabbit  
Isotype:  
IgG

GenBank Accession Number:  
GeneID (NCBI):  
24796  
Full Name:  
sialophorin  
Calculated MW:  
38 kDa

Purification Method:  
Affinity purification  
CloneNo.:  
W3/13  
Excitation/Emission maxima  
wavelengths:  
495 nm / 524 nm

## Applications

Tested Applications:  
FC  
Species Specificity:  
Rat


## Storage

Storage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> rat splenocytes were surface stained with 0.25 ug FITC Plus Anti-Rat CD43 (W3/13) Rabbit Recombinant Antibody (FITC-65556, Clone: W3/13), or FITC-conjugated Rabbit IgG control Rabbit PolyAb (FITC-30000). Cells were not fixed.