

For Research Use Only

FITC Plus Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: FITC-65571



Basic Information

Catalog Number:

FITC-65571

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Protein A purification

Size:

100ug , 500 µg/ml

Full Name:

CloneNo.:

MOPC-173

Source:

Mouse

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Isotype:

IgG2a, kappa

Applications

Tested Applications:

FC

Species Specificity:

n/a

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with 0.2 μ g FITC Plus Mouse IgG2a Isotype Control (MOPC-173) RecAb (FITC-65571, Clone: MOPC-173) (red), or unstained cells (blue), or 5 μ l FITC Plus Anti-Human CD4 (OKT4) Mouse IgG2a RecAb (FITC-65603, Clone: OKT4) (green). Cells were not fixed. Lymphocytes were gated.