

For Research Use Only

# MultiPro™ 5CFLX Anti-Human IL-3RA (Polyclonal)



Catalog Number: G13655-1-5C

## Basic Information

Catalog Number:  
G13655-1-5C

Size:  
10ug, Concentration: 500ug/mL by  
Bradford method using BSA as the  
standard;

Source:  
Rabbit

Isotype:  
IgG

Immunogen Catalog Number:  
AG4452

GenBank Accession Number:  
BC035407

GeneID (NCBI):  
3563

ENSEMBL Gene ID:  
ENSG00000185291

UNIPROT ID:  
P26951

Full Name:  
MultiPro™ 5CFLX Anti-Human IL-3RA  
(Polyclonal)

Conjugate:  
5CFLX

Full Oligo Sequence:  
CGGAGATGTGATAAGAGACAGTTCT  
CTGGCCAGGCTCCCATATAAGAAA

Barcode Sequence:  
TTCTCTGGCCAGGCT

## Applications

Tested Applications:  
Single Cell (Intra)

Species Specificity:  
Human

## Background Information


## Storage

Storage:  
2-8°C  
Storage Buffer:  
PBS with 1mM EDTA and 0.09% sodium azide


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



G13655-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G13655-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.