

For Research Use Only

MultiPro™ 5CFLX Anti-Human Ezrin (Polyclonal)



Catalog Number: G26056-1-5C

Basic Information

Catalog Number:

G26056-1-5C

Size:

10ug, Concentration: 500ug/mL by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23530

GenBank Accession Number:

BC013903

GeneID (NCBI):

7430

ENSEMBL Gene ID:

ENSG00000092820

UNIPROT ID:

P15311

Full Name:

MultiPro™ 5CFLX Anti-Human Ezrin
(Polyclonal)

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGAGACAGGCCA
CCAGTGGCTCCCCATATAAGAAA

Barcode Sequence:

GCCACCAGTGGCTCC

Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

Human

Background Information

Storage

Storage:

2-8°C

Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G26056-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G26056-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.