

chromotek® HA-Peptide

Product Code: ap

Product Information

Description: HPLC-purified HA-peptide

Sequence: YPYDVPDYA

Molecular Weight: 1.1 kDa

Form: HPLC-purified peptide, lyophilized powder

Purity: 95 %

Application: Elution of HA-tagged fusion proteins bound to HA-Traps (ata, atma, atd)

Reconstitution: Reconstitute peptide in 1550 µl PBS, which results in a final concentration of 500 µM. Vortex for 1 min to dissolve the powder.

Shipment: Shipped at ambient temperature

Storage Condition: Upon receipt store at -20°C. Aliquot upon reconstitution and store at -20°C. Avoid repeated freeze/thaw cycles.

Stability: Stable for 1 year upon receipt

Product Sizes

Product	Product Code	Size
HA-Peptide	ap	1 mg

Related Products

Product	Code
HA-Trap Agarose	ata
HA-Trap Magnetic Agarose	atma
HA-Trap Magnetic Particles M-270	atd
HA-Trap Agarose Kit	atak
HA-Trap Magnetic Agarose Kit	atmak
HA-Trap Magnetic Particles M-270 Kit	atdk
ChromoTek HA Monoclonal Antibody (7C9)	7c9
HA Tag Polyclonal Antibody	51064-2-AP
HA Tag Monoclonal Antibody	66006-2-Ig
HA Tag Recombinant Antibody	81290-1-RR

chromotek® HA-Peptide

Product Code: ap

Contact

Proteintech North America (HQ)

Proteintech Group, Inc.
5500 Pearl Street, Suite 400
Rosemont, IL 60018 USA

1-888-472-4522

proteintech@ptglab.com

Proteintech Europe

Transmission (6th Fl)
6 Atherton Street
M3 3GS, Manchester, UK

+44 161 839 3007

europe@ptglab.com

ChromoTek & Proteintech Germany

Fraunhoferstr. 1
82152, Planegg-Martinsreid
Germany

+49 89 124 148 850

germany@ptglab.com

Disclaimer

For research use only (RUO), not for diagnostic or therapeutic use.

ChromoTek® is a registered trademark of ChromoTek GmbH, part of Proteintech Group.