

For Research Use Only

HRP-conjugated ABCE1 Polyclonal antibody

Catalog Number: HRP-28548

Featured Product



Basic Information

Catalog Number:

HRP-28548

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29724

GenBank Accession Number:

BC016283

GeneID (NCBI):

6059

UNIPROT ID:

P61221

Full Name:

ATP-binding cassette, sub-family E (OABP), member 1

Calculated MW:

599 aa, 67 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB

Species Specificity:

Human

Positive Controls:

WB : A549 cells, HEK-293 cells, HeLa cells, K-562 cells, MCF-7 cells

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with HRP-28548 (ABCE1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.