

# IHCeasy<sup>®</sup> Desmocollin 2 Ready-To-Use IHC Kit

Catalog Number: **KHC0091**

## General Information

Sample type:  
FFPE tissue  
Cited sample type:  
Reactivity:  
Human  
Cited Reactivity:

Assay type:  
Immunohistochemistry  
Primary antibody type:  
Rabbit Polyclonal  
Secondary antibody type:  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

## Background

Desmocollin 2 (DSC2) is a member of the desmocollin subfamily of the cadherin superfamily. The desmocollins (DSC1-3), along with the desmogleins (DSG1-4), are found primarily in epithelia and cardiac muscle where they constitute the adhesive proteins of the intercellular desmosome junctions and are required for cell adhesion and desmosome formation. Mutations in DSC2 gene are associated with arrhythmogenic right ventricular dysplasia-11.


## Synonyms

ARVD11, Cadherin family member 2, CDHF2, desmocollin 2, Desmocollin 3, Desmosomal glycoprotein II, Desmosomal glycoprotein III, DG2, DGII/III, DSC2, DSC3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using KHC0091 (Desmocollin 2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using KHC0091 (Desmocollin 2 IHC Kit).