

IHC*easy* ASUN Ready-To-Use IHC Kit

Catalog Number: **KHC0400**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

| Component | Size | Concentration |
|--------------------------|--------------------|---------------|
| Antigen Retrieval Buffer | 100 mL | 50× |
| Washing Buffer | 100 mL ×2 | 20× |
| Blocking Buffer | 5 mL | RTU |
| Primary Antibody | 5 mL | RTU |
| Secondary Antibody | 5 mL | RTU |
| Chromogen Component A | 0.2 mL | RTU |
| Chromogen Component B | 4 mL | RTU |
| Signal Enhancer | 5 mL | RTU |
| Counter Staining Reagent | 5 mL | RTU |
| Mounting Media | 5 mL | RTU |
| Control Slide | 1 slide (Optional) | FFPE |
| Datasheet | 1 Copy | |
| Manual | 1 Copy | |

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

Asunder (ASUN) is a conserved protein that has been reported as a critical regulator of dynein localization during *Drosophila* spermatogenesis. Mammalian and *Drosophila* ASUN proteins revealed that they are 43% identical and 64% similar. Human ASUN (also known as GCT1 or C12ORF11) promotes dynein recruitment and centrosomal tethering to the nucleus at mitotic entry.


Synonyms

ASUN, C12orf11, GCT1, NET48, NY SAR 95, Sarcoma antigen NY SAR 95


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human liver tissue slide using KHC0400 (ASUN IHC Kit).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC0400 (ASUN IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC0400 (ASUN IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC0400 (ASUN IHC Kit).