

IHC*easy* GRHPR Ready-To-Use IHC Kit

Catalog Number: **KHC0520**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

GRHPR (Glyoxylate reductase/hydroxypyruvate reductase) is also named as GLXR and belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. It has glyoxylate reductase (GR), hydroxypyruvate reductase (HPR) and D-glycerate dehydrogenase (DGDH) activities. Defects in GRHPR are the cause of hyperoxaluria primary type 2 (HP2). It can exist as a homodimer.


Synonyms

GLXR, GLYD, GRHPR, PH2


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC0520 (GRHPR IHC Kit).




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC0520 (GRHPR IHC Kit).




Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC0520 (GRHPR IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC0520 (GRHPR IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC0520 (GRHPR IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC0520 (GRHPR IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC0520 (GRHPR IHC Kit).