

IHC*easy* NUFIP1 Ready-To-Use IHC Kit

Catalog Number: **KHC1984**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Fragile X syndrome, the most common cause of inherited mental retardation, is caused by the absence of FMRP (Fragile X Mental Retardation Protein). FMRP is an RNA binding protein reported to be involved in translational control, notably at postsynaptic sites of protein synthesis as a part of a multiprotein/mRNA complex. NUFIP1 is one of the several FMRP-interacting proteins. NUFIP can act as a pol II-specific basal transcriptional activator in vitro and when ectopically overexpressed in vivo. NUFIP can directly activate promoters by enhancing the ATP-dependent release of hyperphosphorylated form of pol II from open transcription complexes.


Synonyms

NUFIP, Nuclear FMRP-interacting protein 1, FMR1-interacting protein NUFIP1


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung tissue slide using KHC1984 (NUFIP1 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse thymus tissue slide using KHC1984 (NUFIP1 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat thymus tissue slide using KHC1984 (NUFIP1 IHC Kit).