

IHC*easy* DMRT1 Ready-To-Use IHC Kit

Catalog Number: **KHC2027**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

DMRT1, also named DMT1, belongs to the DMRT family. DMRT1 is one of the most conserved transcription factors in a variety of species and is involved in sex determination and germ cell differentiation. DMRT1 prevents the sexual switch from male to female, represses pluripotency, and balances mitosis versus meiosis. DMRT1 also maintains spermatogenesis and replenishes mGSCs by collaborating with Plzf and Sal-like protein 4.


Synonyms

DMRT1, DMT1

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC2027 (DMRT1 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC2027 (DMRT1 IHC Kit).