

IHC*easy* ANKRD17 Ready-To-Use IHC Kit

Catalog Number: **KHC2108**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

ANKRD17 belongs to the family of ankyrin repeat-containing proteins, and contains two distinct arrays of ankyrin repeats in its amino-terminal region. It also contains a nuclear export signal, nuclear localization signal, and a cyclin-binding RXL motif. Localization of this protein to the nucleus has been shown experimentally, and interactions between this protein and cyclin-dependent kinase 2 have been observed.


Synonyms

ANKRD17, ankyrin repeat domain 17, GTAR, KIAA0697, NY BR 16


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using KHC2108 (ANKRD17 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse cerebellum tissue slide using KHC2108 (ANKRD17 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC2108 (ANKRD17 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue slide using KHC2108 (ANKRD17 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC2108 (ANKRD17 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC2108 (ANKRD17 IHC Kit).