

# IHC*easy* BIRC7 Ready-To-Use IHC Kit

Catalog Number: **KHC2781**

## General Information

Sample type:  
FFPE tissue  
Cited sample type:  
Reactivity:  
Human, Mouse, Rat  
Cited Reactivity:

Assay type:  
Immunohistochemistry  
Primary antibody type:  
Rabbit Polyclonal  
Secondary antibody type:  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

## Background

BIRC7 also known as Livin, is one of the family members of the inhibitor of apoptosis protein (IAP). BIRC7 contains a BIR domain that has been shown to bind and inhibit caspases-3, -7 and -9, playing a key role in inhibiting apoptosis. Overexpression of BIRC7 in tumors has been linked with increased chemo- and radio-resistance, recurrence, and decreased patient survival. BIRC7 is closely associated with melanoma disease progression, apoptosis resistance, and chemotherapeutic drug resistance.


## Synonyms

Baculoviral IAP repeat-containing protein 7, Baculoviral IAP repeat-containing protein 7 30kDa subunit, EC:2.3.2.27, KIAP, Livin


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2781 (BIRC7 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using KHC2781 (BIRC7 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using KHC2781 (BIRC7 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using KHC2781 (BIRC7 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse skin tissue slide using KHC2781 (BIRC7 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat heart tissue slide using KHC2781 (BIRC7 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC2781 (BIRC7 IHC Kit).