

# IHCeasy<sup>®</sup> PCGF6 Ready-To-Use IHC Kit

Catalog Number: **KHC2886**

## General Information

Sample type:  
FFPE tissue  
Cited sample type:  
Reactivity:  
Human, Mouse, Rat  
Cited Reactivity:

Assay type:  
Immunohistochemistry  
Primary antibody type:  
Rabbit Polyclonal  
Secondary antibody type:  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

## Background

PCGF6, also named as MBLR or RNF134, is a 350 amino acid protein, which contains one RING-type zinc finger. PCGF6 localizes in the nucleus and is widely expressed in many tissues. PCGF6 as a transcriptional repressor may modulate the levels of histone H3K4Me3 by activating KDM5D histone demethylase.


## Synonyms

MBLR, Mel18 and Bmi1-like RING finger, PCGF 6, RNF134


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2886 (PCGF6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC2886 (PCGF6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC2886 (PCGF6 IHC Kit).