

IHC*easy* DDX4/VASA Ready-To-Use IHC Kit

Catalog Number: **KHC2943**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

DEAD box proteins are characterized by nine conserved sequence motifs located on two functional domains. Domain I contains six of these motifs, including the Q motif and the Walker A motif, motifs Ia and Ib, the Walker B motif, and motif III, which may act to link ATPase and helicase activities of the protein. DDX4, a member of the DEAD box family of ATP-dependent RNA helicases, plays a central role in several aspects of germ cell development. Its function is not only required during gametogenesis in the adult but is also essential for the specification of the germ cell lineage during embryogenesis.


Synonyms

DDX4, EC:3.6.4.13, Probable ATP-dependent RNA helicase DDX4, VASA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC2943 (DDX4/VASA IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC2943 (DDX4/VASA IHC Kit).