

IHCeasy[®] CWC15 Ready-To-Use IHC Kit

Catalog Number: **KHC2982**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

CWC15, also named as C11orf5 or Spliceosome-associated protein CWC15 homolog, is a 229 amino acid protein, which belongs to the CWC15 family. CWC15 as a Component of the PRP19-CDC5L complex that forms an integral part of the spliceosome and is required for activating pre-mRNA splicing.


Synonyms

AD-002, C11orf5, Cwf15, HSPC148, ORF5


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2982 (CWC15 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC2982 (CWC15 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using KHC2982 (CWC15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC2982 (CWC15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC2982 (CWC15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat pancreas tissue slide using KHC2982 (CWC15 IHC Kit).