

IHCeasy[®] TAOK2 Ready-To-Use IHC Kit

Catalog Number: **KHC3031**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Recombinant
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

TAOK2(Serine/threonine-protein kinase TAO2) is also named as KIAA0881, MAP3K17, PSK, PSK1 and belongs to the protein kinase superfamily. TAOK2 is known to modulate gene transcription via its activation of several MAP kinase pathways, including c-Jun amino-terminal kinase (JNK), and has been implicated in the modulation of actin and microtubule dynamics. The time course of TAOK2 expression is consistent with a role for it in neuronal differentiation.


Synonyms

EC:2.7.11.1, hKFC-C, KIAA0881, MAP3K17, Prostate-derived STE20-like kinase 1


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC3031 (TAOK2 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC3031 (TAOK2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC3031 (TAOK2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC3031 (TAOK2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC3031 (TAOK2 IHC Kit).