

IHCeasy[®] NSFL1C Ready-To-Use IHC Kit

Catalog Number: **KHC3045**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Synonyms

NSFL1 cofactor p47, p97 cofactor p47, UBX domain-containing protein 2C, UBXN2C


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC3045 (NSFL1C IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC3045 (NSFL1C IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat ovary tissue slide using KHC3045 (NSFL1C IHC Kit).