

IHCeasy[®] UCHL5 Ready-To-Use IHC Kit

Catalog Number: **KHC3078**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

UCHL5 (Ubiquitin carboxyl-terminal hydrolase isozyme L5) is also named as UCH37 and belongs to the peptidase C12 family. It is functionally linked with multiple protein complexes and signal transduction pathways.


Synonyms

UCH37, AD-019, CGI 70, CGI70, CGI-70


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC3078 (UCLH5 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC3078 (UCLH5 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using KHC3078 (UCLH5 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using KHC3078 (UCLH5 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat pancreas tissue slide using KHC3078 (UCLH5 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC3078 (UCLH5 IHC Kit).