

# PDGFR beta Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00006-1

## Capture Antibody Information

Catalog Number:  
82943-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
230214C10  
Reactivity:  
human

Conjugate:  
Unconjugated  
Full name:  
platelet-derived growth factor  
receptor, beta polypeptide  
Gene ID:  
5159

## Detection Antibody Information

Catalog Number:  
82943-4-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
230214E10  
Reactivity:  
human  
GenBank:  
BC032224

Conjugate:  
Unconjugated  
Full name:  
platelet-derived growth factor  
receptor, beta polypeptide  
Gene ID:  
5159

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-10 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP00006-1 targets PDGFR beta in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDGFRB Recombinant antibody, PBS Only (Capture) 82943-2-PBS (230214C10). 100 µg.  
Concentration 1 mg/mL.

Detection antibody: PDGFRB Recombinant antibody, PBS Only (Detector) 82943-4-PBS (230214E10). 100 µg.  
Concentration 1 mg/mL.

Alternative PDGFR beta matched antibody pairs: MP00006-4, MP01260-1, MP01260-2, MP01260-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00006-1, PDGFR beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82943-2-PBS. Detection antibody: 82943-4-PBS. Standard: Eg0147. Range: 0.313-10 ng/mL.