

# EIF2S1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00076-2

## Capture Antibody Information

**Catalog Number:**  
82936-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230245D1  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag1645

**Conjugate:**  
Unconjugated  
**Full name:**  
eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa  
**Gene ID:**  
1965

## Detection Antibody Information

**Catalog Number:**  
82936-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230245E7  
**Reactivity:**  
human  
**GenBank:**  
BC002513  
**Immunogen Catalog Number:**  
Ag1645

**Conjugate:**  
Unconjugated  
**Full name:**  
eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa  
**Gene ID:**  
1965

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00076-2 targets EIF2S1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EIF2S1 Recombinant antibody, PBS Only (Capture) 82936-4-PBS (230245D1). 100 µg. Concentration 1 mg/mL.

Detection antibody: EIF2S1 Recombinant antibody, PBS Only (Detector) 82936-5-PBS (230245E7). 100 µg. Concentration 1 mg/mL.

Alternative EIF2S1 matched antibody pairs: MP00076-1, MP00076-3, MP00076-4, MP50179-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00076-2, EIF2S1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82936-4-PBS. Detection antibody: 82936-5-PBS. Standard: Ag1645. Range: 0.625-20 ng/mL