

For Research Use Only



ABCA13 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00084-2

Capture Antibody Information

Catalog Number:
83060-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230288H4
Reactivity:
Human

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 13
Gene ID:
154664

Detection Antibody Information

Catalog Number:
83060-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230288G4
Reactivity:
Human
GenBank:
NM_152701

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 13
Gene ID:
154664

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00084-2 targets ABCA13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABCA13 Recombinant antibody, PBS Only (Capture) 83060-3-PBS (230288H4). 100 µg. Concentration 1 mg/ml.

Detection antibody: ABCA13 Recombinant antibody, PBS Only (Detector) 83060-2-PBS (230288G4). 100 µg. Concentration 1 mg/ml.

Alternative ABCA13 matched antibody pairs: MP00084-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00084-2, ABCA13
Recombinant Matched Antibody Pair - BSA and
Azide free, measured by cytometric bead array.
Capture antibody: 83060-3-PBS. Detection
antibody: 83060-2-PBS. Standard: SY00509. Range:
0.312-20 ng/mL