

For Research Use Only



# TNFR1/CD120a Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00099-2

## Capture Antibody Information

Catalog Number:  
83080-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
230322D3  
Reactivity:  
human

Conjugate:  
Unconjugated  
Full name:  
tumor necrosis factor receptor  
superfamily, member 1A  
Gene ID:  
7132

## Detection Antibody Information

Catalog Number:  
83080-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
230322E8  
Reactivity:  
human  
GenBank:  
BC010140

Conjugate:  
Unconjugated  
Full name:  
tumor necrosis factor receptor  
superfamily, member 1A  
Gene ID:  
7132

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00099-2 targets TNFR1/CD120a in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFR1/CD120a Recombinant antibody, PBS Only (Capture) 83080-1-PBS (230322D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: TNFR1/CD120a Recombinant antibody, PBS Only (Detector) 83080-2-PBS (230322E8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00099-2, TNFR1/CD120a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83080-1-PBS. Detection antibody: 83080-2-PBS. Standard: Eg0531. Range: 1.25-20 ng/mL