

For Research Use Only



SLC35E4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00107-1

Capture Antibody Information

Catalog Number:
82999-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240089C12

Reactivity:
Human

Immunogen Catalog Number:
Ag34360

Conjugate:
Unconjugated

Full name:
solute carrier family 35, member E4

Gene ID:
339665

Detection Antibody Information

Catalog Number:
82999-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240089A12

Reactivity:
Human

GenBank:
BC040191

Immunogen Catalog Number:
Ag34360

Conjugate:
Unconjugated

Full name:
solute carrier family 35, member E4

Gene ID:
339665

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00107-1 targets SLC35E4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC35E4 Recombinant antibody, PBS Only (Capture) 82999-2-PBS (240089C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC35E4 Recombinant antibody, PBS Only (Detector) 82999-1-PBS (240089A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00107-1, SLC35E4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82999-2-PBS. Detection antibody: 82999-1-PBS. Standard: Ag34360. Range: 1.25-20 ng/mL