

uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00129-1

Capture Antibody Information

Catalog Number:

82953-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230148B10

Reactivity:

human, mouse

Immunogen Catalog Number:

Ag28056

Conjugate:

Unconjugated

Full name:

PLAUR

Gene ID:

5329

Detection Antibody Information

Catalog Number:

82953-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230148E10

Reactivity:

human

GenBank:

BC002788

Immunogen Catalog Number:

Ag28056

Conjugate:

Unconjugated

Full name:

PLAUR

Gene ID:

5329

Applications

Tested Applications:

Sandwich ELISA

Range:

125-8000 pg/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00129-1 targets uPAR/CD87 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: uPAR/CD87 Recombinant antibody, PBS Only (Capture/Detector) 82953-2-PBS (230148B10). 100 µg. Concentration 1 mg/mL.

Detection antibody: uPAR/CD87 Recombinant antibody, PBS Only (Detector) 82953-3-PBS (230148E10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00129-1, uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82953-2-PBS. Detection antibody: 82953-3-PBS. Standard: Ag28056. Range: 125-8000 pg/mL