

For Research Use Only

# SALL4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00137-1

## Capture Antibody Information

Catalog Number:

83039-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230284D4

Reactivity:

Human

Immunogen Catalog Number:

Ag16076

Conjugate:

Unconjugated

Full name:

sal-like 4 (Drosophila)

Gene ID:

57167

## Detection Antibody Information

Catalog Number:

83039-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230284B7

Reactivity:

Human

GenBank:

BC111714

Immunogen Catalog Number:

Ag16076

Conjugate:

Unconjugated

Full name:

sal-like 4 (Drosophila)

Gene ID:

57167

## Applications

Tested Applications:

Sandwich ELISA

Range:

0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00137-1 targets SALL4 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: SALL4 Recombinant antibody, PBS Only (Capture) 83039-2-PBS (230284D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: SALL4 Recombinant antibody, PBS Only (Detector) 83039-1-PBS (230284B7). 100 µg. Concentration 1 mg/mL.

Alternative SALL4 matched antibody pairs: MP00137-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00137-1, SALL4 Recombinant Matched Antibody Pair - BSA and Azide free, measured by sandwich ELISA. Capture antibody: 83039-2-PBS. Detection antibody: 83039-1-PBS. Standard: Ag16076. Range: 0.156-10 ng/mL