

For Research Use Only

SLC25A28 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00142-2

Capture Antibody Information

Catalog Number:

83127-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230392C1

Reactivity:

Human

Immunogen Catalog Number:

Ag34507

Conjugate:

Unconjugated

Full name:

solute carrier family 25, member 28

Gene ID:

81894

Detection Antibody Information

Catalog Number:

83127-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230392F8

Reactivity:

Human

GenBank:

BC076399

Immunogen Catalog Number:

Ag34507

Conjugate:

Unconjugated

Full name:

solute carrier family 25, member 28

Gene ID:

81894

Applications

Tested Applications:

Cytometric bead array

Range:

1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00142-2 targets SLC25A28 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A28 Recombinant antibody, PBS Only (Capture/Detector) 83127-2-PBS (230392C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC25A28 Recombinant antibody, PBS Only (Detector) 83127-5-PBS (230392F8). 100 µg. Concentration 1 mg/ml.

Alternative SLC25A28 matched antibody pairs: MP00142-1, MP00142-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00142-2, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-2-PBS. Detection antibody: 83127-5-PBS. Standard: Ag34507. Range: 1.25-20 ng/mL