

For Research Use Only

SLC25A28 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00142-3

Capture Antibody Information

Catalog Number:
83127-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230392A3

Reactivity:
Human

Immunogen Catalog Number:
Ag34507

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 28

Gene ID:
81894

Detection Antibody Information

Catalog Number:
83127-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230392C5

Reactivity:
Human

GenBank:
BC076399

Immunogen Catalog Number:
Ag34507

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 28

Gene ID:
81894

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00142-3 targets SLC25A28 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A28 Recombinant antibody, PBS Only (Capture) 83127-1-PBS (230392A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC25A28 Recombinant antibody, PBS Only (Detector) 83127-3-PBS (230392C5). 100 µg. Concentration 1 mg/ml.

Alternative SLC25A28 matched antibody pairs: MP00142-1, MP00142-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00142-3, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-1-PBS. Detection antibody: 83127-3-PBS. Standard: Ag34507. Range: 0.625-20 ng/mL.