

CELSR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00165-2**

Capture Antibody Information

Catalog Number:
83115-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230319C10
Reactivity:
human

Conjugate:
Unconjugated
Full name:
cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)
Gene ID:
1952

Detection Antibody Information

Catalog Number:
83115-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230319D10
Reactivity:
human
GenBank:
NM_001408

Conjugate:
Unconjugated
Full name:
cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)
Gene ID:
1952

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00165-2 targets CELSR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: CELSR2 Recombinant antibody, PBS Only (Capture) 83115-2-PBS (230319C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CELSR2 Recombinant antibody, PBS Only (Detector) 83115-3-PBS (230319D10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

