

# GBX1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00201-3**

## Capture Antibody Information

**Catalog Number:**  
83232-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240108H1  
**Reactivity:**  
Human

**Conjugate:**  
Unconjugated  
**Full name:**  
gastrulation brain homeobox 1  
**Gene ID:**  
2636

## Detection Antibody Information

**Catalog Number:**  
83232-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240108A8  
**Reactivity:**  
Human  
**GenBank:**  
NM\_001098834

**Conjugate:**  
Unconjugated  
**Full name:**  
gastrulation brain homeobox 1  
**Gene ID:**  
2636

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00201-3 targets GBX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GBX1 Recombinant antibody, PBS Only (Capture/Detector) 83232-2-PBS (240108H1). 100 µg. Concentration 1 mg/mL.

Detection antibody: GBX1 Recombinant antibody, PBS Only (Detector) 83232-1-PBS (240108A8). 100 µg. Concentration 1 mg/mL.

Alternative GBX1 matched antibody pairs: MP00201-1, MP00201-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00201-3, GBX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83232-2-PBS. Detection antibody: 83232-1-PBS. Standard: SY01055. Range: 1.25-20 ng/mL.