

For Research Use Only



MKNK2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00216-4

Capture Antibody Information

Catalog Number:
83235-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240045E9
Reactivity:
human
Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated
Full name:
MAP kinase interacting
serine/threonine kinase 2
Gene ID:
2872

Detection Antibody Information

Catalog Number:
83235-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240045D11
Reactivity:
human
GenBank:
BC018345
Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated
Full name:
MAP kinase interacting
serine/threonine kinase 2
Gene ID:
2872

Applications

Tested Applications:
Sandwich ELISA

Range:
0.313-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00216-4 targets MKNK2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: MKNK2 Recombinant antibody, PBS Only (Capture) 83235-5-PBS (240045E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: MKNK2 Recombinant antibody, PBS Only (Capture/Detector) 83235-4-PBS (240045D11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00216-4, Human MKNK2 Recombinant Matched Antibody Pair - PBS only. 83235-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10546. 83235-4-PBS was HRP conjugated as the detection antibody. Range: 0.313-10 ng/mL.



The mean MKNK2 concentration was determined to be 13.14 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load.