

For Research Use Only



IL-17RB Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00240-1

Capture Antibody Information

Catalog Number:
83149-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230439H7

Reactivity:
human

Conjugate:
Unconjugated

Full name:
interleukin 17 receptor B

Gene ID:
55540

Detection Antibody Information

Catalog Number:
83149-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230439E1

Reactivity:
human

GenBank:
NM_018725

Conjugate:
Unconjugated

Full name:
interleukin 17 receptor B

Gene ID:
55540

Applications

Tested Applications:
Cytometric bead array

Range:
0.078-2.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00240-1 targets IL-17RB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-17RB Recombinant antibody, PBS Only (Capture) 83149-2-PBS (230439H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-17RB Recombinant antibody, PBS Only (Detector) 83149-1-PBS (230439E1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00240-1, IL-17RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83149-2-PBS. Detection antibody: 83149-1-PBS. Standard: SY00888. Range: 0.078-2.5 ng/mL