

B7-H3/CD276 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00248-1**

Capture Antibody Information

Catalog Number:
83152-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230305F11
Reactivity:
human

Conjugate:
Unconjugated
Full name:
CD276 molecule
Gene ID:
80381

Detection Antibody Information

Catalog Number:
83152-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230305A5
Reactivity:
human
GenBank:
NM_001024736

Conjugate:
Unconjugated
Full name:
CD276 molecule
Gene ID:
80381

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00248-1 targets B7-H3/CD276 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD276 Recombinant antibody, PBS Only (Capture) 83152-4-PBS (230305F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD276 Recombinant antibody, PBS Only (Detector) 83152-1-PBS (230305A5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00248-1, B7-H3/CD276 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83152-4-PBS. Detection antibody: 83152-1-PBS. Standard: Eg0075. Range: 0.625-20 ng/mL.