

N4BP3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00260-1

Capture Antibody Information

Catalog Number:

83178-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230547E2

Reactivity:

Human

Immunogen Catalog Number:

Ag10194

Conjugate:

Unconjugated

Full name:

Nedd4 binding protein 3

Gene ID:

23138

Detection Antibody Information

Catalog Number:

83178-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230547C6

Reactivity:

Human

GenBank:

BC053323

Immunogen Catalog Number:

Ag10194

Conjugate:

Unconjugated

Full name:

Nedd4 binding protein 3

Gene ID:

23138

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00260-1 targets N4BP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: N4BP3 Recombinant antibody, PBS Only (Capture) 83178-4-PBS (230547E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: N4BP3 Recombinant antibody, PBS Only (Capture/Detector) 83178-3-PBS (230547C6). 100 µg. Concentration 1 mg/mL.

Alternative N4BP3 matched antibody pairs: MP00260-2, MP00260-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00260-1, N4BP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83178-4-PBS. Detection antibody: 83178-3-PBS. Standard: Ag10194. Range: 0.625-20 ng/mL