

# CEACAM8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00293-3**

## Capture Antibody Information

**Catalog Number:**  
83105-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230353A10  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
carcinoembryonic antigen-related cell adhesion molecule 8  
**Gene ID:**  
1088

## Detection Antibody Information

**Catalog Number:**  
83105-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230353B7  
**Reactivity:**  
human  
**GenBank:**  
BC026263

**Conjugate:**  
Unconjugated  
**Full name:**  
carcinoembryonic antigen-related cell adhesion molecule 8  
**Gene ID:**  
1088

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00293-3 targets CEACAM8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CEACAM8 Recombinant antibody, PBS Only (Capture) 83105-2-PBS (230353A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CEACAM8 Recombinant antibody, PBS Only (Detector) 83105-1-PBS (230353B7). 100 µg. Concentration 1 mg/mL.

Alternative CEACAM8 matched antibody pairs: MP00293-1, MP00293-2, MP00293-4, MP50039-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00293-3, CEACAM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83105-2-PBS. Detection antibody: 83105-1-PBS. Standard: Eg32133. Range: 0.625-20 ng/mL.