

CHD6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00298-2

Capture Antibody Information

Catalog Number:
83122-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230397D6
Reactivity:
Human

Conjugate:
Unconjugated
Full name:
chromodomain helicase DNA binding protein 6
Gene ID:
84181

Detection Antibody Information

Catalog Number:
83122-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230397B12
Reactivity:
Human
GenBank:
NM_032221

Conjugate:
Unconjugated
Full name:
chromodomain helicase DNA binding protein 6
Gene ID:
84181

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00298-2 targets CHD6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHD6 Recombinant antibody, PBS Only (Capture) 83122-6-PBS (230397D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: CHD6 Recombinant antibody, PBS Only (Detector) 83122-4-PBS (230397B12). 100 µg. Concentration 1 mg/mL.

Alternative CHD6 matched antibody pairs: MP00298-1, MP00298-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00298-2, CHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83122-6-PBS. Detection antibody: 83122-4-PBS. Standard: SY01008. Range: 0.156-5 ng/mL.