

For Research Use Only

CA15-3,MUC1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00335-2

Capture Antibody Information

Catalog Number:
83311-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240190F1

Reactivity:
Human

Immunogen Catalog Number:
Ag20366

Conjugate:
Unconjugated

Full name:
mucin 1, cell surface associated

Gene ID:
4582

Detection Antibody Information

Catalog Number:
83311-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240190A2

Reactivity:
Human

GenBank:
BC120975

Immunogen Catalog Number:
Ag20366

Conjugate:
Unconjugated

Full name:
mucin 1, cell surface associated

Gene ID:
4582

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00335-2 targets CA15-3,MUC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CA15-3,MUC1 Recombinant antibody, PBS Only (Capture) 83311-2-PBS (240190F1). 100 µg. Concentration 1 mg/ml.

Detection antibody: CA15-3,MUC1 Recombinant antibody, PBS Only (Detector) 83311-1-PBS (240190A2). 100 µg. Concentration 1 mg/ml.

Alternative CA15-3,MUC1 matched antibody pairs: MP00335-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00335-2, CA15-3, MUC1 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83311-2-PBS. Detection antibody: 83311-1-PBS. Standard: Ag20366. Range: 1.56-100 ng/mL.