

SAT1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00344-3**

Capture Antibody Information

Catalog Number:
83319-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240157H3
Reactivity:
human
Immunogen Catalog Number:
Ag1154

Conjugate:
Unconjugated
Full name:
spermidine/spermine N1-acetyltransferase 1
Gene ID:
6303

Detection Antibody Information

Catalog Number:
83319-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240157G2
Reactivity:
human
GenBank:
BC008424
Immunogen Catalog Number:
Ag1154

Conjugate:
Unconjugated
Full name:
spermidine/spermine N1-acetyltransferase 1
Gene ID:
6303

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00344-3 targets SAT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SAT1 Recombinant antibody, PBS Only (Capture/Detector) 83319-2-PBS (240157H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: SAT1 Recombinant antibody, PBS Only (Capture/Detector) 83319-1-PBS (240157G2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00344-3, SAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83319-2-PBS. Detection antibody: 83319-1-PBS. Standard:null. Range: null.