

# CMTM6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00352-1**

## Capture Antibody Information

**Catalog Number:**  
83334-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240194B3  
**Reactivity:**  
Human  
**Immunogen Catalog Number:**  
Ag34445

**Conjugate:**  
Unconjugated  
**Full name:**  
CKLF-like MARVEL transmembrane domain containing 6  
**Gene ID:**  
54918

## Detection Antibody Information

**Catalog Number:**  
83334-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240194A12  
**Reactivity:**  
Human  
**GenBank:**  
BC002797  
**Immunogen Catalog Number:**  
Ag34445

**Conjugate:**  
Unconjugated  
**Full name:**  
CKLF-like MARVEL transmembrane domain containing 6  
**Gene ID:**  
54918

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
2.5-80 ng/mL (Cytometric Bead Array)  
**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00352-1 targets CMTM6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CMTM6 Recombinant antibody, PBS Only (Capture) 83334-2-PBS (240194B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: CMTM6 Recombinant antibody, PBS Only (Detector) 83334-1-PBS (240194A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00352-1, CMTM6 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83334-2-PBS. Detection antibody: 83334-1-PBS. Standard: Ag34445. Range: 2.5-80 ng/mL.