

BAP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00357-1**

Capture Antibody Information

Catalog Number:
83371-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240349A7
Reactivity:
human
Immunogen Catalog Number:
Ag34415

Conjugate:
Unconjugated
Full name:
BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)
Gene ID:
8314

Detection Antibody Information

Catalog Number:
83371-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240349F5
Reactivity:
human
GenBank:
BC001596
Immunogen Catalog Number:
Ag34415

Conjugate:
Unconjugated
Full name:
BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)
Gene ID:
8314

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00357-1 targets BAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BAP1 Recombinant antibody, PBS Only (Capture) 83371-1-PBS (240349A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: BAP1 Recombinant antibody, PBS Only (Detector) 83371-2-PBS (240349F5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00357-1, BAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83371-1-PBS. Detection antibody: 83371-2-PBS. Standard: Ag34415. Range: 1.25-80 ng/mL.