

CNKS2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00389-2

Capture Antibody Information

Catalog Number:
83354-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240350C8
Reactivity:
Human
Immunogen Catalog Number:
Ag21174

Conjugate:
Unconjugated
Full name:
connector enhancer of kinase suppressor of Ras 2
Gene ID:
22866

Detection Antibody Information

Catalog Number:
83354-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240350H8
Reactivity:
Human
GenBank:
BC126121
Immunogen Catalog Number:
Ag21174

Conjugate:
Unconjugated
Full name:
connector enhancer of kinase suppressor of Ras 2
Gene ID:
22866

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00389-2 targets CNKS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CNKS2 Recombinant antibody, PBS Only (Capture) 83354-1-PBS (240350C8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CNKS2 Recombinant antibody, PBS Only (Detector) 83354-4-PBS (240350H8). 100 µg. Concentration 1 mg/mL.

Alternative CNKS2 matched antibody pairs: MP00389-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00389-2, CNKSR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83354-1-PBS. Detection antibody: 83354-4-PBS. Standard: Ag21174. Range: 6.25-100ng/mL.