

For Research Use Only

# Rat MCP-1/CCL2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00390-1

## Capture Antibody Information

Catalog Number:  
83370-3-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240341E10  
Reactivity:  
rat

Conjugate:  
Unconjugated  
Full name:  
chemokine (C-C motif) ligand 2  
Gene ID:  
24770

## Detection Antibody Information

Catalog Number:  
83370-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240341D11  
Reactivity:  
rat  
GenBank:  
NM\_031530

Conjugate:  
Unconjugated  
Full name:  
chemokine (C-C motif) ligand 2  
Gene ID:  
24770

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00390-1 targets MCP-1/CCL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat MCP-1/CCL2 Recombinant antibody, PBS Only (Capture) 83370-3-PBS (240341E10). 100 µg. Concentration 1 mg/mL.

Detection antibody: Rat MCP-1/CCL2 Recombinant antibody, PBS Only (Detector) 83370-1-PBS (240341D11). 100 µg. Concentration 1 mg/mL.

Alternative MCP-1/CCL2 matched antibody pairs: MP00168-1, MP00168-2, MP00168-3, MP00168-4, MP00390-2, MP00390-3, MP00938-1, MP00938-2, MP00938-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00390-1, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-1-PBS. Standard: Eg0397. Range: 0.313-20 ng/mL