

SLC6A7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00399-1**

Capture Antibody Information

Catalog Number:

83412-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240395D5

Reactivity:

human

Immunogen Catalog Number:

Ag34413

Conjugate:

Unconjugated

Full name:

solute carrier family 6 (neurotransmitter transporter, L-proline), member 7

Gene ID:

6534

Detection Antibody Information

Catalog Number:

83412-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240395E7

Reactivity:

human

GenBank:

BC093785

Immunogen Catalog Number:

Ag34413

Conjugate:

Unconjugated

Full name:

solute carrier family 6 (neurotransmitter transporter, L-proline), member 7

Gene ID:

6534

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00399-1 targets SLC6A7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC6A7 Recombinant antibody, PBS Only (Capture) 83412-3-PBS (240395D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC6A7 Recombinant antibody, PBS Only (Detector) 83412-2-PBS (240395E7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00399-1, SLC6A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83412-3-PBS. Detection antibody: 83412-2-PBS. Standard: Ag34413. Range: 0.625-80 ng/mL