

# SLC25A17 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00409-1**

## Capture Antibody Information

**Catalog Number:**  
83404-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240382G8

**Reactivity:**  
Human

**Immunogen Catalog Number:**  
Ag34310

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 25 (mitochondrial carrier; peroxisomal membrane protein, 34kDa), member 17

**Gene ID:**  
10478

## Detection Antibody Information

**Catalog Number:**  
83404-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240382E8

**Reactivity:**  
Human

**GenBank:**  
BC012998

**Immunogen Catalog Number:**  
Ag34310

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 25 (mitochondrial carrier; peroxisomal membrane protein, 34kDa), member 17

**Gene ID:**  
10478

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
2.5-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00409-1 targets SLC25A17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A17 Recombinant antibody, PBS Only (Capture/Detector) 83404-2-PBS (240382G8). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC25A17 Recombinant antibody, PBS Only (Detector) 83404-1-PBS (240382E8). 100 µg. Concentration 1 mg/mL.

Alternative SLC25A17 matched antibody pairs: MP00409-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00409-1, SLC25A17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83404-2-PBS. Detection antibody: 83404-1-PBS. Standard: Ag34310. Range: 2.5-80 ng/mL.